



ERBAdry[®]

The NEW generation of dry solvents



DASITGROUP

CARLO ERBA


REAGENTS

ERBAdry[®] : NEW generation of caps to prevent from moisture and to improve customer experience

- **100ml & 1l** amber glass bottles to protect chemicals from UV. Some solvents are also available on **molecular sieve**
- **FLIP TOP CAP** always ready to be closed
- **2 temper evident rings** showing the double unused product warranty
- **3 large areas** available to improve usability
- **Septum double layers** to prevent chemical from contamination
- **2 types of septa** selected for better chemical compatibility :
 - Orange septum : Butyl rubber / PTFE
 - Black septum : Bromobutyl / PTFE



Evolution of water content



Water content for Chloroform stabilized with amylene, performed with Karl Fischer method over 6 weeks. Up to 28 samplings from the same bottle using a 18 Gauge needle.

To find out more 



Find out the list of solvents in our ERBAdry® range.

Produit	Application	Code 100 ml	Code 1 l	Code 1l with sieve
Acetone	RS ERBAdry®	P00510D03	P00510D16	
Acetonitrile	RS ERBAdry®	P00610D03	P00610D16	P00610TD16
Butanol-1	RS ERBAdry®	P01710D03	P01710D16	
tert-Butanol	RS ERBAdry®	P01910D03	P01910D16	
Butyle acetate	RS ERBAdry®	P00110D03	P00110D16	
tert-Butylmethylether	RS ERBAdry®	P09210D03	P09210D16	
Chloroform	RS ERBAdry® stabilized with amylene	P02410AD03	P02410AD16	
Chloroform	RS ERBAdry® stabilized with ethanol	P02410ED03	P02410ED16	
Cyclohexane	RS ERBAdry®	P02510D03	P02510D16	
Dichloromethane	RS ERBAdry® stabilized with amylene	P02910AD03	P02910AD16	P02910ATD16
Dichloromethane	RS ERBAdry® stabilized with ethanol	P02910ED03	P02910ED16	
Diethyl ether	RS ERBAdry® stabilized with BHT	P04410D03	P04410D16	P04410TD16
1,2-Dimethoxyethane	RS ERBAdry®			P03010TD16
N,N - Dimethylformamide	RS ERBAdry®	P03410D03	P03410D16	P03410TD16
Dimethylsulfoxide	RS ERBAdry®	5081D03	5081D16	P03510TD16
1,4-Dioxane	RS ERBAdry® stabilized with BHT	P03610D03	P03610D16	
Ethanol absolu	RS ERBAdry®	P013A10D03	P013A10D16	
Ethyle acetate	RS ERBAdry®	P00210D03	P00210D16	P00210TD16
Formamide	RS ERBAdry®	P61510D03	P61510D16	
n-Heptane	RS ERBAdry®	P05010D03	P05010D16	
n-Hexane 99%	RS ERBAdry®	P05230D03	P05230D16	
n-Hexane 95%	RS ERBAdry®	P05210D03	P05210D16	
Isopentane	RS ERBAdry®	P06510D03	P06510D16	
Isopropyl ether	RS ERBAdry®	P04310D03	P04310D16	
Methanol	RS ERBAdry®	P09310D03	P09310D16	
Methylcyclohexane	RS ERBAdry®	P05810D03	P05810D16	
Methylethylketone	RS ERBAdry®	P02010D03	P02010D16	
Methylisobutylketone	RS ERBAdry®	P06010D03	P06010D16	
N-Methylpyrrolidone-2	RS ERBAdry®	P08710D03	P08710D16	
4-Methyltetrahydropyran	RS ERBAdry®	P99910D03	P99910D16	
n-Pentane	RS ERBAdry®	P06410D03	P06410D16	
Propanol-1	RS ERBAdry®	P09410D03	P09410D16	
Propanol-2	RS ERBAdry®	P09510D03	P09510D16	
Pyridine	RS ERBAdry®	P06710D03	P06710D16	
Tetrahydrofuran	RS ERBAdry® stabilized with BHT	P07010D03	P07010D16	P07010TD16
Toluene	RS ERBAdry®	P07110D03	P07110D16	P07110TD16

Check the availability on our website
www.carloerbareagents.com





ITALIA

CARLO ERBA Reagents S.r.l.
Via Raffaele Merendi 22
20007 Cornaredo (MI)

Servizio Clienti

servizioclienticer@dgroup.it
Tel.: +39 02 93 99 190
Fax: +39 02 93 99 10 01



FRANCE

CARLO ERBA Reagents SAS
Chaussée du Vexin,
Parc d'affaire des Portes
27106 Val de Reuil

Service Client

serviceclient@cer.dgroup.it
Tél.: +33 2 32 09 20 00
Fax: +33 2 32 59 11 89



DEUTSCHLAND

CARLO ERBA Reagents GmbH
Denzlinger Str. 27
79312 Emmendingen

Kundendienst

info.de@cer.dgroup.it
Tel.: +49 07641 46 881 90
Fax: +49 07641 46 881 919



ESPAÑA

CARLO ERBA Reagents S.A.
Calle Filadors 35
6ª Planta 5ª Puerta
08208 Sabadell (BCN)

Servicio Cliente

serviciocliente@cer.dgroup.it
Tel.: +34 93 693 37 35
Fax: +34 93 724 31 68



ALL OTHER COUNTRIES

Customer Service
export@cer.dgroup.it
Tel.: +33 2 32 09 20 00
Fax: +33 2 32 59 11 89

